

# Product Release MiniBath™

Part Number: 120342


Considerable improvements were made to the MiniBath module of the  $\mu$ DISS ProfilerPLUS™ including \*8\* individual sample stations and computer controlled stirring.

## $\mu$ DISS ProfilerPLUS Product Features:

- Real time concentration monitoring
- Low sample consumption (< 0.5 mg)
- Media sparing
- No filtration required
- 6 or 8 UV fiber optic probes
- Temperature controlled
- Applications
  - Solubility
  - Dissolution
  - Polymorph Screening
  - Salt Selection
  - Excipient measurements



## Product Improvements:



- Small volume dissolution made easy
- 8 sample stations provided
- Each station has individual stirring control
- Computer controlled stirring speed
- Staggered starts for each station
- Countdown for introduction of media solution

## Ordering information:

Part Number: 120342

Customers with a six channel MiniBath may trade in their MiniBath for the eight channel MiniBath. MiniBath sold separately or as a component of the  $\mu$ DISS ProfilerPLUS

